

Statistical Studies: Experiment

Luke A. Siemens*
(Dated: October 3, 2016)

XXX

I. INTRODUCTION

The first component of this experiment focus on the statistical atool sand teqniques for analyzing adata and statistical distribution. Spesificaly methods for fitting statistical distributions and teqniques for determining if data is consistent with a theoretical model are investigated. In the second component of theis experiment the statistical method investigated in the first section are applied to the practical problem of calibrating labritory equipment, and error analysis of the subsequent mesurements. In the corse of this lab we will be investigating the Poisson distribution, the Gaussian distributheon the χ^2 distribution and the χ^2 statistic and theoir application to the analysis of counting statistics and the callibration of a multi-channel analyzer.

II. THEORY

The statistic of counting random events is modeled by assuming the process in question can be described as a large number of simple binary events which each have a small probability of leading to a count being mesured, where the number of binary event n is $n \sim 1/x$ each event having a probability of the binary even ocuring is $p \sim 1/x$ for $x \ll 1$. The probability that ν events will occur in a given interval is described by the binomial distribution,

$$B_{n,p}(\nu) = \binom{n}{\nu} p^\nu (1-p)^{n-\nu}$$

In the limiting case of a continuous process $p \rightarrow 0$ and $n \rightarrow \infty$ the distribution simplifies to the Poisson distribution.

$$P_\mu(\nu) = \frac{e^{-\mu} \mu^\nu}{\nu!} \quad (1)$$

where the expected number of events $\bar{\nu}$ is $\bar{\nu} = \sum_{\nu=0}^{\infty} \nu P_\mu(\nu) = \mu$ and the variance is $\sigma_\nu^2 = \sum_{\nu=0}^{\infty} (\nu - \mu)^2 P_\mu(\nu) = \mu$. The Poisson distribution can be simplified in the limiting case where μ is large, in that case the Poisson distribution (1) aproches a gaussian distribution with mean μ and variance μ . In this limit the distribution becomes,

$$P_\mu(\nu) \approx G_{\mu, \sqrt{\mu}}(\nu) = \frac{1}{\sqrt{2\pi\mu}} e^{-(\nu-\mu)^2/2\mu} \quad (2)$$

A. χ^2 Test

The χ^2 test test is a method for determining the probability that a set of mesurements is consistent with a given model assuming Gaussian errors. The χ^2 test uses the χ^2 statistic which in the case of descreat variables is defined as,

$$\chi^2 = \sum_i \left(\frac{O_i - E_i}{\sigma_i} \right)^2$$

where O_i is the observed value, E_i is the expected value and σ_i is the standard deviation. Given the expected value is distributed as a Poisson distribution with mean $\mu = E_i$ and varianece $\sigma_i = \sqrt{E_i}$ then the χ^2 statistic is,

$$\chi^2 = \sum_i \frac{(O_i - E_i)^2}{E_i} \quad (3)$$

Assuming the errors are gaussaing then the χ^2 statistic is described by the associated χ^2 distribution. The χ^2 distribution which has the cumulative distribution function (CDF),

$$\text{Prob}_d(\chi^2 \leq \chi_o^2) = \frac{1}{2^{d/2} \Gamma(d/2)} \int_0^{\chi_o^2} t^{d/2-1} e^{-t/2} dt \quad (4)$$

So the probability of measureing a χ^2 as large or larger than χ_o^2 is $\text{Prob}_d(\chi^2 \geq \chi_o^2) = 1 - \text{Prob}_d(\chi^2 \leq \chi_o^2)$.

B. Poisson errors and modified χ^2 statistic

The statistic given by equation (3) can be used with the χ^2 assuming O_i is described by a gaussian distribution but that assumption is invalid when E_i is not sufficently large that the Poisson distribution can be aproximated as a Gaussian distribution. If when E_i is small the χ^2 statistic must be modiefied to account for the asymmetry of the Poisson distribution. For the gaussian distribution the probability of getting a value below the lower bound of a one sigma confidence interval is $p = \frac{1}{2} - \text{erf}(\frac{1}{\sqrt{2}})/2$,

* lsiemens@uvic.ca

for the modified χ^2 statistic instead of the standard deviation lets use the value $|\nu_o - \mu|$ where ν_o solves the equation $\text{Prob}_\mu(\nu \leq \nu_o) = p$ or $\text{Prob}_\mu(\nu \geq \nu_o) = 1 - p$ depending on whether the observed value is above or below the expected value. To evaluate these probabilities for μ . The CDF of the poisson distribution is,

$$\text{Prob}_\mu(\nu \leq \nu_o) = \sum_{k=0}^{\nu_o} \frac{e^{-\mu} \mu^k}{k!}$$

where ν_o is a positive integer. When $\mu \sim 1$ a continuous analogue to the Poisson CDF is needed to evaluate ν_o for the adjusted χ^2 statistic. In the paper *Continuous Counterparts of Poisson and Binomial Distributions and their Properties*[1] A. Ilien defines a continuous analogue to the Poisson distribution, the CDF of this distribution is

$$\text{Prob}_\mu(\nu \leq \nu_o) = \frac{\Gamma(\nu_o + 1, \mu)}{\Gamma(\nu_o + 1)}$$

where ν_o is a positive real number, $\Gamma(x, \lambda)$ is the incomplete Gamma function $\Gamma(x, \lambda) = \int_\lambda^\infty t^{x-1} e^{-t} dt$ and $\Gamma(x) = \Gamma(x, 0)$ is the Gamma function. Note that if ν_o is a positive integer then,

$$\text{Prob}_\mu(\nu \leq \nu_o) = \frac{\Gamma(\nu_o + 1, \mu)}{\Gamma(\nu_o + 1)} = \sum_{k=0}^{\nu_o} \frac{e^{-\mu} \mu^k}{k!}$$

Using the continuous analogue to the Poisson distribution the adjusted χ^2 statistic is defined as,

$$\chi^2 = \sum_i \left(\frac{O_i - E_i}{\nu_{o_i} - E_i} \right)^2$$

$$\frac{\Gamma(\nu_{o_i} + 1, E_i)}{\Gamma(\nu_{o_i} + 1)} = \begin{cases} p & \text{if } O_i < E_i \\ 1 - p & \text{if } O_i > E_i \end{cases}$$

where $p = \frac{1}{2} - \text{erf}(\frac{1}{\sqrt{2}})/2$. It should be noted that for $E_i > 30$ the relative difference between $|\nu_{o_i} - E_i|$ and $\sqrt{E_i}$ is less than ten percent and that if $E_i < 0.08503$ then the lower ν_{o_i} will be such that $\nu_{o_i} > E_i$.

C. Energy resolution

The resolution R of a spectrometer is,

$$R = \frac{\delta E}{E} \cdot 100 \quad (5)$$

where δE is the Full Width at Half Maximum (FWHM) of an energy peak and where E is the energy at the maximum of the energy peak. The resolution R gauges the detectors ability to resolve spectral features.

III. DESIGN AND RESULTS

The experiment is split into two components, in the first section the distribution of counts produced by a sintilator(s) is measured and compared to theoretical distributions, in the second section the methods investigated for analyzing distributions is used to calibrate a multi-channel analyzer and to measure the resolution of the device.

A. Section I

The distribution produced by one sintilator for the cases of $\mu \approx 5$ and $\mu \approx 100$ where measured. We setup the sintilator with a $-1.5kV$ bias voltage. The signal from the sintilator was passed into an amplifier and then through a discriminator. the signal from the discriminator was then passed to a timer/counter. Two sets of counts were collected with this setup each with 500 measurements. In the first set the discriminator cutoff was tuned to produce a low count rate $\mu \approx 5$, in the second set the cutoff was tuned to produce a high count rate of $\mu \approx 100$.

In the second part two sintilators were used to reject any signal that was not correlated between the devices, so as to reduce the random noise from the sintilators. The sintilators were stacked vertically with approximately 6cm vertical separation. The signals from the two sintilators were then passed into two amplifiers and then into two discriminators. The two signals were then passed into a correlator with correlation setting set to two. The output from the correlator was then connected to a counter/timer. As in the first part two sets of data was taken with 500 measurements each, one with $\mu \approx 5$ and the other with $\mu \approx 100$.

B. Section II

using a NaI detector we calibrated a multi-channel analyzer by fitting the energy peak of known isotopes, using the analysis techniques investigated in Section I. Using the calibrated measurements of a known sample to find the energy resolution of the multi-channel analyzer.

IV. ANALYSIS

V. CONCLUSION

-
- [1] A. Ilienکو, *Continuous Counterparts of Poisson and Binomial Distributions and their Properties*, (Annales Univ. Sci. Budapest., Sect. Comp. 39(2013) 137-147).